





NLRP3 Antibody

Product Code	CSB-PA015871LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P20
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 3 protein (30-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 3 (Angiotensin/vasopressin receptor AII/AVP-like) (Caterpillar protein 1.1) (CLR1.1) (Cold-induced autoinflammatory syndrome 1 protein) (Cryopyrin) (PYRIN-containing APAF1-like protein 1), NLRP3, C1orf7 CIAS1 NALP3 PYPAF1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRP3

Image



Immunohistochemistry of paraffin-embedded human bladder cancer using CSB-PA015871LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA015871LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)