



# NLRP3 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA015871LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96P20
<b>Immunogen</b>	Recombinant Human NACHT, LRR and PYD domains-containing protein 3 protein (30-102AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	NACHT, LRR and PYD domains-containing protein 3 (Angiotensin/vasopressin receptor All/AVP-like) (Caterpillar protein 1.1) (CLR1.1) (Cold-induced autoinflammatory syndrome 1 protein) (Cryopyrin) (PYRIN-containing APAF1-like protein 1), NLRP3, C1orf7 CIAS1 NALP3 PYPAF1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Immunology
<b>Target Names</b>	NLRP3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.