



PDGFRA Antibody, Biotin conjugated

Product Code	CSB-PA017712LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16234
Immunogen	Recombinant Human Platelet-derived growth factor receptor alpha protein (30-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet-derived growth factor receptor alpha (PDGF-R-alpha) (PDGFR-alpha) (EC 2.7.10.1) (Alpha platelet-derived growth factor receptor) (Alpha-type platelet-derived growth factor receptor) (CD140 antigen-like family member A) (CD140a antigen) (Platelet-derived growth factor alpha receptor) (Platelet-derived growth factor receptor 2) (PDGFR-2) (CD antigen CD140a), PDGFRA, PDGFR2 RHEPDGFRA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PDGFRA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.