



PEMT Antibody, Biotin conjugated

Product Code	CSB-PA017780LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBM1
Immunogen	Recombinant Human Phosphatidylethanolamine N-methyltransferase protein (2-49AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylethanolamine N-methyltransferase (PEAMT) (PEMT) (EC 2.1.1.17) (EC 2.1.1.71) (PEMT2) (Phospholipid methyltransferase) (PLMT), PEMT, PEMPT PNMT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PEMT