



PFKM Antibody, Biotin conjugated

Product Code	CSB-PA017822LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08237
Immunogen	Recombinant Human ATP-dependent 6-phosphofructokinase, muscle type protein (166-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent 6-phosphofructokinase, muscle type (ATP-PFK) (PFK-M) (EC 2.7.1.11) (6-phosphofructokinase type A) (Phosphofructo-1-kinase isozyme A) (PFK-A) (Phosphohexokinase), PFKM, PFKX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PFKM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.