





PIGP Antibody

Product Code	CSB-PA017979LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57054
Immunogen	Recombinant Human Phosphatidylinositol N-acetylglucosaminyltransferase subunit P protein (1-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol N-acetylglucosaminyltransferase subunit P (EC 2.4.1.198) (Down syndrome critical region protein 5) (Down syndrome critical region protein C) (Phosphatidylinositol-glycan biosynthesis class P protein) (PIG-P), PIGP, DCRC DSCR5 DSCRC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PIGP

Image



Immunohistochemistry of paraffin-embedded human stomach tissue using CSB-PA017979LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA017979LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)