



PIGP Antibody, FITC conjugated

Product Code	CSB-PA017979LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57054
Immunogen	Recombinant Human Phosphatidylinositol N-acetylglucosaminyltransferase subunit P protein (1-282AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol N-acetylglucosaminyltransferase subunit P (EC 2.4.1.198) (Down syndrome critical region protein 5) (Down syndrome critical region protein C) (Phosphatidylinositol-glycan biosynthesis class P protein) (PIG-P), PIGP, DCRC DSCR5 DSCRC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PIGP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.