



PRDX3 Antibody, Biotin conjugated

Product Code	CSB-PA018656LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30048
Immunogen	Recombinant Human Thioredoxin-dependent peroxide reductase, mitochondrial protein (63-256AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin-dependent peroxide reductase, mitochondrial (EC 1.11.1.15) (Antioxidant protein 1) (AOP-1) (HBC189) (Peroxiredoxin III) (Prx-III) (Peroxiredoxin-3) (Protein MER5 homolog), PRDX3, AOP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PRDX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.