



RAB7A Antibody, Biotin conjugated

Product Code	CSB-PA019219LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51149
Immunogen	Recombinant Human Ras-related protein Rab-7a protein (110-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-7a, RAB7A, RAB7
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	RAB7A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.