



SCGB1A1 Antibody, Biotin conjugated

Product Code	CSB-PA020811LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11684
Immunogen	Recombinant Human Uteroglobin protein (22-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uteroglobin (Clara cell phospholipid-binding protein) (CCPBP) (Clara cells 10 kDa secretory protein) (CC10) (Secretoglobin family 1A member 1) (Urinary protein 1) (UP-1) (UP1) (Urine protein 1), SCGB1A1, CC10 CCSP UGB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SCGB1A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.