



# SCGB1A1 Antibody, Biotin conjugated

Product Code	CSB-PA020811LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11684
Immunogen	Recombinant Human Uteroglobulin protein (22-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Uteroglobulin (Clara cell phospholipid-binding protein) (CCPBP) (Clara cells 10 kDa secretory protein) (CC10) (Secretoglobulin family 1A member 1) (Urinary protein 1) (UP-1) (UP1) (Urine protein 1), SCGB1A1, CC10 CCSP UGB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SCGB1A1