



SGK1 Antibody, FITC conjugated

Product Code	CSB-PA021189LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00141
Immunogen	Recombinant Human Serine/threonine-protein kinase Sgk1 protein (247-406AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Sgk1 (EC 2.7.11.1) (Serum/glucocorticoid-regulated kinase 1), SGK1, SGK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SGK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.