



SH2D1A Antibody, HRP conjugated

Product Code	CSB-PA021209LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60880
Immunogen	Recombinant Human SH2 domain-containing protein 1A protein (11-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH2 domain-containing protein 1A (Duncan disease SH2-protein) (Signaling lymphocytic activation molecule-associated protein) (SLAM-associated protein) (T-cell signal transduction molecule SAP), SH2D1A, DSHP SAP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SH2D1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.