



SH2D1A Antibody, HRP conjugated

Product Code	CSB-PA021209LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60880
Immunogen	Recombinant Human SH2 domain-containing protein 1A protein (11-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH2 domain-containing protein 1A (Duncan disease SH2-protein) (Signaling lymphocytic activation molecule-associated protein) (SLAM-associated protein) (T-cell signal transduction molecule SAP), SH2D1A, DSHP SAP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SH2D1A