



SH2D1B Antibody, HRP conjugated

Product Code	CSB-PA021210LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14796
Immunogen	Recombinant Human SH2 domain-containing protein 1B protein (1-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH2 domain-containing protein 1B (EWS/FLI1-activated transcript 2) (EAT-2), SH2D1B, EAT2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SH2D1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.