



SH3BP5 Antibody, HRP conjugated

Product Code	CSB-PA021226LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60239
Immunogen	Recombinant Human SH3 domain-binding protein 5 protein (49-455AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-binding protein 5 (SH3BP-5) (SH3 domain-binding protein that preferentially associates with BTK), SH3BP5, SAB
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SH3BP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.