



SH3BP5 Antibody, FITC conjugated

Product Code	CSB-PA021226LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60239
Immunogen	Recombinant Human SH3 domain-binding protein 5 protein (49-455AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-binding protein 5 (SH3BP-5) (SH3 domain-binding protein that preferentially associates with BTK), SH3BP5, SAB
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SH3BP5