



SLC18A1 Antibody, FITC conjugated

Product Code	CSB-PA021426LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54219
Immunogen	Recombinant Human Chromaffin granule amine transporter protein (66-120AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chromaffin granule amine transporter (Solute carrier family 18 member 1) (Vesicular amine transporter 1) (VAT1), SLC18A1, VAT1 VMAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC18A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.