



# SLC35A2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA021584LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P78381
<b>Immunogen</b>	Recombinant Human UDP-galactose translocator protein (136-234AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	UDP-galactose translocator (Solute carrier family 35 member A2) (UDP-galactose transporter) (UDP-Gal-Tr) (UGT), SLC35A2, UGALT UGT UGTL
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	SLC35A2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.