



SLC35A2 Antibody, Biotin conjugated

Product Code	CSB-PA021584LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78381
Immunogen	Recombinant Human UDP-galactose translocator protein (136-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-galactose translocator (Solute carrier family 35 member A2) (UDP-galactose transporter) (UDP-Gal-Tr) (UGT), SLC35A2, UGALT UGT UGTL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC35A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.