



SLC37A1 Antibody, FITC conjugated

Product Code	CSB-PA021610LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57057
Immunogen	Recombinant Human Glucose-6-phosphate exchanger SLC37A1 protein (245-302AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-6-phosphate exchanger SLC37A1 (Glycerol-3-phosphate permease) (G-3-P permease) (Solute carrier family 37 member 1), SLC37A1, G3PP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC37A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.