



SLC37A1 Antibody, Biotin conjugated

Product Code	CSB-PA021610LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57057
Immunogen	Recombinant Human Glucose-6-phosphate exchanger SLC37A1 protein (245-302AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-6-phosphate exchanger SLC37A1 (Glycerol-3-phosphate permease) (G-3-P permease) (Solute carrier family 37 member 1), SLC37A1, G3PP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC37A1