



# SLC9A1 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA021726LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P19634
<b>Immunogen</b>	Recombinant Human Sodium/hydrogen exchanger 1 protein (36-102AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Sodium/hydrogen exchanger 1 (APNH) (Na(+)/H(+)) antiporter, amiloride-sensitive) (Na(+)/H(+)) exchanger 1) (NHE-1) (Solute carrier family 9 member 1), SLC9A1, APNH1 NHE1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SLC9A1