



SLC9A1 Antibody, HRP conjugated

Product Code	CSB-PA021726LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19634
Immunogen	Recombinant Human Sodium/hydrogen exchanger 1 protein (36-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/hydrogen exchanger 1 (APNH) (Na(+)/H(+)) antiporter, amiloride-sensitive) (Na(+)/H(+)) exchanger 1) (NHE-1) (Solute carrier family 9 member 1), SLC9A1, APNH1 NHE1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC9A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.