




SLFN11 Antibody

Product Code	CSB-PA021760LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7L1
Immunogen	Recombinant Human Schlafen family member 11 protein (670-821AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:100-1:500, IF:1:50-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Schlafen family member 11, SLFN11
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLFN11

Image



Western Blot

Positive WB detected in: Mouse skeletal muscle tissue


All lanes: SLFN11 antibody at 2.8μg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 103 kDa

Observed band size: 103 kDa



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA021760LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA021760LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)