



SMS Antibody, HRP conjugated

Product Code	CSB-PA021854LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52788
Immunogen	Recombinant Human Spermine synthase protein (16-191AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermine synthase (SPMSY) (EC 2.5.1.22) (Spermidine aminopropyltransferase), SMS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SMS