





SNAI2 Antibody

Product Code	CSB-PA021869LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43623
Immunogen	Recombinant Human Zinc finger protein SNAI2 protein (12-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein SNAI2 (Neural crest transcription factor Slug) (Protein snail homolog 2), SNAI2, SLUG SLUGH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNAI2

Image



Immunohistochemistry of paraffin-embedded human gastric tissue using CSB-PA021869LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA021869LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.