





SPINK1 Antibody

Product Code	CSB-PA022575LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00995
Immunogen	Recombinant Human Serine protease inhibitor Kazal-type 1 protein (19-79AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor Kazal-type 1 (Pancreatic secretory trypsin inhibitor) (Tumor-associated trypsin inhibitor) (TATI), SPINK1, PSTI
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPINK1


Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA022575LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA022575LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA022575LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.