



SPINK1 Antibody, HRP conjugated

Product Code	CSB-PA022575LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00995
Immunogen	Recombinant Human Serine protease inhibitor Kazal-type 1 protein (19-79AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor Kazal-type 1 (Pancreatic secretory trypsin inhibitor) (Tumor-associated trypsin inhibitor) (TATI), SPINK1, PSTI
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPINK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.