



STMN1 Antibody, HRP conjugated

Product Code	CSB-PA022857LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16949
Immunogen	Recombinant Human Stathmin protein (13-78AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stathmin (Leukemia-associated phosphoprotein p18) (Metablastin) (Oncoprotein 18) (Op18) (Phosphoprotein p19) (pp19) (Prosolin) (Protein Pr22) (pp17), STMN1, C1orf215 LAP18 OP18
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	STMN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.