



TGFBR1 Antibody, Biotin conjugated

Product Code	CSB-PA023451LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36897
Immunogen	Recombinant Human TGF-beta receptor type-1 protein (157-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TGF-beta receptor type-1 (TGFR-1) (EC 2.7.11.30) (Activin A receptor type II-like protein kinase of 53kD) (Activin receptor-like kinase 5) (ALK-5) (ALK5) (Serine/threonine-protein kinase receptor R4) (SKR4) (TGF-beta type I receptor) (Transforming growth factor-beta receptor type I) (TGF-beta receptor type I) (TbetaR-I), TGFBR1, ALK5 SKR4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TGFBR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.