



TLL1 Antibody, HRP conjugated

Product Code	CSB-PA023595LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43897
Immunogen	Recombinant Human Tolloid-like protein 1 protein (522-643AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tolloid-like protein 1 (EC 3.4.24.-), TLL1, TLL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TLL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.