



TLR3 Antibody, Biotin conjugated

Product Code	CSB-PA023602LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15455
Immunogen	Recombinant Human Toll-like receptor 3 protein (135-297AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Toll-like receptor 3 (CD antigen CD283), TLR3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TLR3