





TOR1A Antibody

Product Code	CSB-PA024067LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14656
Immunogen	Recombinant Human Torsin-1A protein (88-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Torsin-1A (Dystonia 1 protein) (Torsin ATPase-1A) (EC 3.6.4.-) (Torsin family 1 member A), TOR1A, DQ2 DYT1 TA TORA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TOR1A

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA024067LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA024067LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.