



TOR1A Antibody, Biotin conjugated

Product Code	CSB-PA024067LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14656
Immunogen	Recombinant Human Torsin-1A protein (88-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Torsin-1A (Dystonia 1 protein) (Torsin ATPase-1A) (EC 3.6.4.-) (Torsin family 1 member A), TOR1A, DQ2 DYT1 TA TORA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TOR1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.