



TTLL8 Antibody, HRP conjugated

Product Code	CSB-PA025265LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6PVC2
Immunogen	Recombinant Human Protein monoglycylase TTLL8 protein (671-822AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein monoglycylase TTLL8 (EC 6.3.2.-) (Tubulin--tyrosine ligase-like protein 8), TTLL8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTLL8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.