



UBA6 Antibody, HRP conjugated

Product Code	CSB-PA025420LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0AVT1
Immunogen	Recombinant Human Ubiquitin-like modifier-activating enzyme 6 protein (1-44AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like modifier-activating enzyme 6 (Ubiquitin-activating enzyme 6) (EC 6.2.1.45) (Monocyte protein 4) (MOP-4) (Ubiquitin-activating enzyme E1-like protein 2) (E1-L2), UBA6, MOP4 UBE1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.