



ACTG1 Antibody, HRP conjugated

Product Code	CSB-PA02619B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63261
Immunogen	Recombinant Human Actin, cytoplasmic 2 protein (13-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin, cytoplasmic 2 (Gamma-actin) [Cleaved into: Actin, cytoplasmic 2, N-terminally processed], ACTG1, ACTG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.