



YES1 Antibody, HRP conjugated

Product Code	CSB-PA026253LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07947
Immunogen	Recombinant Human Tyrosine-protein kinase Yes protein (1-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase Yes (EC 2.7.10.2) (Proto-oncogene c-Yes) (p61-Yes), YES1, YES
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	YES1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.