





ZMPSTE24 Antibody

Product Code	CSB-PA026504LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75844
Immunogen	Recombinant Human CAAX prenyl protease 1 homolog protein (217-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CAAX prenyl protease 1 homolog (EC 3.4.24.84) (Farnesylated proteins-converting enzyme 1) (FACE-1) (Prenyl protein-specific endoprotease 1) (Zinc metalloproteinase Ste24 homolog), ZMPSTE24, FACE1 STE24
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ZMPSTE24


Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA026504LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA026504LA01HU at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA026504LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.