



ZMPSTE24 Antibody, Biotin conjugated

Product Code	CSB-PA026504LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75844
Immunogen	Recombinant Human CAAX prenyl protease 1 homolog protein (217-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CAAX prenyl protease 1 homolog (EC 3.4.24.84) (Farnesylated proteins-converting enzyme 1) (FACE-1) (Prenyl protein-specific endoprotease 1) (Zinc metalloproteinase Ste24 homolog), ZMPSTE24, FACE1 STE24
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ZMPSTE24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.