



MCL1 Antibody, HRP conjugated

Product Code	CSB-PA03829B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07820
Immunogen	Recombinant Human Induced myeloid leukemia cell differentiation protein Mcl-1 protein (19-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Induced myeloid leukemia cell differentiation protein Mcl-1 (Bcl-2-like protein 3) (Bcl2-L-3) (Bcl-2-related protein EAT/mcl1) (mcl1/EAT), MCL1, BCL2L3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MCL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.