



# MCL1 Antibody, Biotin conjugated

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA03829D0Rb   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| Uniprot No.         | Q07820  |
| Immunogen           | Recombinant Human Induced myeloid leukemia cell differentiation protein Mcl-1 protein (19-244AA)  |
| Raised In           | Rabbit  |
| Species Reactivity  | Human   |
| Tested Applications | ELISA   |
| Form                | Liquid  |
| Conjugate           | Biotin  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| Purification Method | >95%, Protein G purified  |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Alias               | Induced myeloid leukemia cell differentiation protein Mcl-1 (Bcl-2-like protein 3) (Bcl2-L-3) (Bcl-2-related protein EAT/mcl1) (mcl1/EAT), MCL1, BCL2L3 |
| Immunogen Species   | Homo sapiens (Human)  |
| Research Area       | Cell Biology  |
| Target Names        | MCL1  |