



FCGR2A Antibody, HRP conjugated

Product Code	CSB-PA04459B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12318
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor II-a protein (241-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor II-a (IgG Fc receptor II-a) (CDw32) (Fc-gamma RII-a) (Fc-gamma-RIIa) (FcRII-a) (CD antigen CD32), FCGR2A, CD32 FCG2 FCGR2A1 IGFR2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCGR2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.