



LASP1 Antibody, Biotin conjugated

Product Code	CSB-PA04629D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14847
Immunogen	Recombinant Human LIM and SH3 domain protein 1 protein (128-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM and SH3 domain protein 1 (LASP-1) (Metastatic lymph node gene 50 protein) (MLN 50), LASP1, MLN50
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LASP1