



LASP1 Antibody, Biotin conjugated

Product Code	CSB-PA04629D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14847
Immunogen	Recombinant Human LIM and SH3 domain protein 1 protein (128-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM and SH3 domain protein 1 (LASP-1) (Metastatic lymph node gene 50 protein) (MLN 50), LASP1, MLN50
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LASP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.