



CXCL11 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA06119B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O14625 |
| Immunogen | Recombinant Human C-X-C motif chemokine 11 protein (22-94AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | C-X-C motif chemokine 11 (Beta-R1) (H174) (Interferon gamma-inducible protein 9) (IP-9) (Interferon-inducible T-cell alpha chemoattractant) (I-TAC) (Small-inducible cytokine B11), CXCL11, ITAC SCYB11 SCYB9B |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | CXCL11 |