



CXCL11 Antibody, FITC conjugated

Product Code	CSB-PA06119C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14625
Immunogen	Recombinant Human C-X-C motif chemokine 11 protein (22-94AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 11 (Beta-R1) (H174) (Interferon gamma-inducible protein 9) (IP-9) (Interferon-inducible T-cell alpha chemoattractant) (I-TAC) (Small-inducible cytokine B11), CXCL11, ITAC SCYB11 SCYB9B
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CXCL11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.