



FLT1 Antibody, FITC conjugated

Product Code	CSB-PA07269C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17948
Immunogen	Recombinant Human Vascular endothelial growth factor receptor 1 protein (1048-1328AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor receptor 1 (VEGFR-1) (EC 2.7.10.1) (Fms-like tyrosine kinase 1) (FLT-1) (Tyrosine-protein kinase FRT) (Tyrosine-protein kinase receptor FLT) (FLT) (Vascular permeability factor receptor), FLT1, FLT FRT VEGFR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FLT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.