



CD70 Antibody, HRP conjugated

Product Code	CSB-PA08009B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32970
Immunogen	Recombinant Human CD70 antigen protein (62-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD70 antigen (CD27 ligand) (CD27-L) (Tumor necrosis factor ligand superfamily member 7) (CD antigen CD70), CD70, CD27L CD27LG TNFSF7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD70